

Sound SmartDR High Frequency Cesium Technique Chart

Weight	0 - 20 lb	21 - 45 lb	46 - 80 lb	81 - 100 lb	100+ lb	Avian
Abdomen	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/4.0 mAs	80 kVp/7.5 mAs	90 kVp/7.5 mAs	55 kVp/10 mAs
Thorax	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/4.0 mAs	80 kVp/7.5 mAs	90 kVp/7.5 mAs	
Extremity	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	
Pelvis	80 kVp/3.0 mAs	80 kVp/4.0 mAs	80 kVp/7.5 mAs	90 kVp/7.5 mAs	100 kVp/7.5 mAs	
Skull	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	80 kVp/3.0 mAs	
Shoulder	80 kVp/4.0 mAs	80 kVp/4.0 mAs	80 kVp/4.0 mAs	80 kVp/6.0 mAs	80 kVp/6.0 mAs	
Spine	80 kVp/7.5 mAs	80 kVp/7.5 mAs	80 kVp/18.0 mAs	80 kVp/18.0 mAs	80 kVp/30.0 mAs	

Note: If images are "grainy" (under-exposed) INCREASE mAs

If images are "saturated" (over-exposed) DECREASE mAs

*Due to variability in X-Ray generators optimal techniques may vary from the values shown above

YOUR IMAGE IS OUR BUSINESS.

 A VCA COMPANY
 800.819.5538 soundvet.com
 install.soundvet.com



Sound SmartDR High Frequency Cesium Technique Chart

Weight	0 - 20 lb	21 - 45 lb	46 - 80 lb	81 - 100 lb	100+ lb	Avian
Abdomen						
Thorax						
Extremity						
Pelvis						
Skull						
Shoulder						
Spine						

Note: If images are "grainy" (under-exposed) INCREASE mAs

If images are "saturated" (over-exposed) DECREASE mAs

*Due to variability in X-Ray generators optimal techniques may vary from the values shown above

YOUR IMAGE IS OUR BUSINESS.

800.819.5538 soundvet.com
install.soundvet.com

