

eFilm Setup Guide

This guide will walk you through setting up eFilm on your hospitals computers. Below is information you will need to enter after the software has installed. If at any point you require assistance, please contact Technical Support at 1-800-819-5538

AE Title:


IP Address:

Port: 104


IC Port: 3333

1. Double click the **set-up** icon to run the InstallShield Wizard.

NOTE: You must be logged in as an Administrator. In Windows 7/Vista - right click the **setup** icon and choose to "Run as Administrator"



2. Click **Install**. The setup files will extract. This can take several minutes.



3. Click **Next**




4. Enter
- Username: SoundEklin
 - Company Name: SoundEklin



5. Click **Next** to continue the installation.




6. Click Next at Choose Destination Location window




7. Enter the following information
- AE Title: Use the name of the computer you are on (EXAM2, PHARM, TREATMENT, etc)
 - Port: 4006
 - click Next



8. Make sure **Typical** is chosen and click **Next**




9. Click **Next** to start copying files



10. Click **OK** when the InstallShield Wizard tells you that it will need to restart the machine




11. Choose **Yes, I want to restart my computer now** and click **Finish**



12. When your computer restarts, double click the eFilm icon on the desktop




13. Click **Register** when the software starts




14. Click **Yes** when it asks if you want to continue the Registration



15. Choose **Client** when it asks for the Registration Type



16. Enter the IP address of your server in the **Primary Server** field. This is found at the top of this document



17. Click **Verify**. The Status field should read *Passed*

| | Address: | Port: | Status: | |
|-------------------|---------------|-------|---------|---------------------------------------|
| Primary Server: | 192.168.1.201 | 8080 | Passed | <input type="button" value="Verify"/> |
| Secondary Server: | | 8080 | | <input type="button" value="Verify"/> |

18. Click **OK**, then click **Yes** when it tells you the Secondary Server is offline




19. Click **OK** to close the **Client License** window

20. The eFilm software will now load. Click the **Settings** button



21. Ensure the following are checked/unchecked, then click **OK**
- a. **Checked:** Bypass user authentication
 - b. **Checked:** Bypass login when automated
 - c. **Unchecked:** Enable audit
 - d. **Checked:** Disable login timeout
 - e. **Checked:** Overlay Fixed
 - f. **Unchecked:** Limit Features



The image shows a screenshot of the 'Administrative Settings' dialog box. The dialog has a blue title bar with the text 'Administrative Settings' and a close button (X) in the top right corner. The main content area is divided into several sections:


- Login Gatekeepers:** A table with columns 'Description', 'Domain Name / Web Service', and 'Type'. Below the table are input fields for 'Description', 'Domain Name', and 'Type' (with a dropdown menu currently set to 'DOMAIN'). There is a 'Browse...' button next to the 'Domain Name' field. Below these fields are buttons for 'Add', 'Update', 'Delete', and 'Verify...'. There is also a blank space for a file path.
- Login Overrides:** Contains two checked checkboxes: 'Bypass user authentication' and 'Bypass login when automated'. There is also an unchecked checkbox for 'Enable audit' and a 'Channel:' field with the value '8085'.
- Login Timeout:** Contains a text field 'Lock eFilm after the user has been inactive for' followed by a spinner box set to '10' and the word 'minutes'. There is a checked checkbox for 'Disable login timeout'.
- Overlay Setting:** Contains a checked checkbox for 'Overlay Fixed'.
- Limit Features:** Contains an unchecked checkbox for 'Limit eFilm Features' and five checked checkboxes: 'Disable Study Manager', 'Disable Help (via F1)', 'Disable Right-click menu', 'Disable Close Study', and 'Disable Close eFilm'.

An 'OK' button is located in the bottom right corner of the dialog.

22. When the software opens, click **Edit**, then click **Properties**



23. Choose the **Image Channel** tab from the Edit Properties page



24. Enter the following information in the Device section
- a. Description: Fusion Server
 - b. AE Title: See Page 1 (This must be entered exactly as it appears, including case)
 - c. IP Address: See Page 1 for the Server's IP Address
 - d. DICOM Port: 104
 - e. IC Port: 3333
 - f. Timeout: 5
 - g. **Check** the Default box

| Device | | | | | | |
|---------------|--------------|---------------|------------|---------|----------------|---|
| Description | AE Title | IP Address | DICOM Port | IC Port | Timeout (secs) | |
| Fusion Server | VARCE0611C34 | 192.168.0.201 | 104 | 3333 | 5 | <input checked="" type="checkbox"/> Default |

Verify

25. Click the **Add** button

26. Single click the **Fusion Server** in the Servers section

27. Click the **Verify** button. If you do not receive a *verified connection* message, please contact Tech Support at 1-800-819-5538

| | |
|--------|-------------------------------------|
| Verify | Fusion Server - connection verified |
|--------|-------------------------------------|

28. Click the **OK** button to return to the eFilm Study Manager

29. Search for your study using Image Channel in the Study Manager

