

Sound SmartPACS™ - Toolbar Quick Reference Guide



Screen Layout Group

Changes the layout of the images in the viewport. The layout can be set to a maximum of 6x6



Screen Mode

Collapses the thumbnail panel allowing the image to be displayed in the entire viewport



Fill

Fills empty windows in the viewport with images from the study being viewed



Image Scope Group

Indicates which images the manipulation tools will be applied to

Tools are applied to the selected image only

Tools are applied to all images in the selected stack

Tools are applied to all images within all stacks




Overlay

Shows or hides the overlay information for all images in the viewport



Previous/Next Stack

Navigate to the previous or next stack



Rotate 90 Left/Right

Rotates the selected image 90 degrees to the left or right



Vertical/Horizontal Mirror

Flips the selected image on its vertical or horizontal axis



Pan

Moves the selected image around the screen



Window Level Group

Changes the Window and Level of the selected image



Changes the window (brightness) and level (contrast) of the selected image



Applies an optimized window and level setting for the selected image



Enable or disable the Look Up Table (LUT)



Saves the current window and level setting as a window/level preset



Applies a window/level preset to the selected image



Zoom Group

Modifies the magnification of the selected image



Zoom - Zooms the entire selected image



Magnifying Glass - Magnifies a portion of the selected image



Zoom True Size - Displays the selected image with a size ratio of 1:1



Invert

Inverts the Black and White values in the selected image



Revert

Reverts any Pan, Window/Level, and Zoom modifications made to the selected image

Sound SmartPACS™ - Toolbar Quick Reference Guide



Measurement Tool Group - Contains all measurement tools

Linear Measurement Tool - Measures the length of a straight line

Arrow - Draws an arrow on the image

Angle - Measures the angle between two lines that intersect

Cobb Angle - Measures the angle between two lines that do not intersect

Vertebral Heart Score - Custom Measurement that calculates the Vertebral Heart Score

Norberg-Winkel - Custom Measurement that calculates Hip Displasia

Circle - Measures the radius, area, and average intensity of the drawn region

Ellipse - Measures the area and average intensity of the drawn region

3 Point Circle - Draws a circle based on 3 points

Area - Measures the area of a drawn region

Text - Adds a text annotation to the image

Toggles all annotations on/off

Edit a single annotation

Delete a single annotation

Delete all annotations



Export Study
Downloads the selected study to the local disk as a DICOMDIR



Screen Capture
Creates a JPEG file of the selected image including annotations



Viewer Settings
Access the settings for the Image Viewer



Help
Access help files and tutorial videos for SmartPACS

