

Body Part Examined	View Position	SID (in)	SID (cm)	kVp	mAs	mA	Sec
NAVICULAR	D65PR_PADIO (65 Degree)	22"	56cm	80	1.00	20	0.05
	PAPR_PADIO (Skyline)	22"	56cm	80	1.00	20	0.05
	LM	26"	66cm	80	1.00	20	0.05
	DLPMO/DMPLO	26"	66cm	80	1.00	20	0.05
	DP	26"	66cm	80	1.00	20	0.05
FOOT	DLPMO/DMPLO	26"	66cm	80	1.00	20	0.05
	Lm	26"	66cm	80	1.00	20	0.05
	Dp	26"	66cm	80	1.00	20	0.05
	D65PR_PADIO (65 Degree)	22"	56cm	80	1.00	20	0.05
	D65PR_PADIO (65 DegreeLat Obl)	22"	56cm	80	1.00	20	0.05
	D65PR_PADILO (65 Degree Med Obl)	22"	56cm	80	1.00	20	0.05
PASTERN	LM	26"	66cm	80	1.00	20	0.05
	DLPMO/DMPLO	26"	66cm	80	1.00	20	0.05
	Dp	26"	66cm	80	1.00	20	0.05
FETLOCK	LM	26"	66 cm	80	1.00	20	0.05
	DLPMO/DMPLO	26"	66 cm	80	1.00	20	0.05
	DP	26"	66 cm	80	1.00	20	0.05
	Flexed LM	26"	66 cm	80	1.00	20	0.05
	Flexed DP	26"	66 cm	80	1.00	20	0.05
CARPUS	LM	26"	66cm	80	1.00	20	0.05
	DLPMO/DMPLO	26"	66cm	80	1.00	20	0.05
	DP	26"	66cm	80	1.00	20	0.05
	Flexed LM	26"	66cm	80	1.00	20	0.05
	DPR_DDIO (Skyline)	22"	56cm	80	1.00	20	0.05
TARSUS	LM	26"	66cm	80	1.00	20	0.05
	DLPMO/DMPLO	26"	66cm	80	1.00	20	0.05
	DP	26"	66cm	80	1.00	20	0.05
	Flexed LM	26"	66cm	80	1.00	20	0.05
	D_PL (Skyline)	22"	56cm	80	1.00	20	0.05
STIFLE	LM	26"	66cm	80	2.00	20	0.13
	CDL_CRMO	26"	66cm	80	2.00	20	0.13
	CDM_CRLO	26"	66cm	80	2.00	20	0.13
	CD_CR	26"	66cm	80	5.00	10	0.50
	Flexed LM	26"	66cm	80	2.00	20	0.13
SHOULDER	LM	26"	66cm	80	5.00	10	0.50
	CRM_CDLO	26"	66cm	80	5.00	10	0.50
	ML	26"	66cm	80	5.00	10	0.50
CERVICAL SPINE	LAT (Cranial)	26"	66cm	80	2.50	20	0.16
	LAT (Middle)	26"	66cm	80	3.20	20	0.21
	LAT (Caudal)	26"	66cm	80	5.00	10	0.50
SKULL	LAT	26"	66cm	80	2.50	20	0.16
	LED_RTVO	26"	66cm	80	2.50	20	0.16
	RTD_LEVO	26"	66cm	80	2.50	20	0.16
	DV	26"	66cm	80	5.00	10	0.50
WITHERS	LAT	26"	66cm	80	2.50	20	0.50
THORACIC SPINE	LAT	26"	66cm	80	5.00	10	0.50
LUMBAR SPINE	LAT	26"	66cm	80	5.00	10	0.50

