

Agfa CR 15-X Technique Chart - High Frequency

Abdomen

CM	2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
kVp	68	72	76	80	84	88	92	96	100	100	100	100	100	100	100
mAs	7.5	7.5	7.5	7.5	7.5	7.5	7.5	10	10	10	10	12.5	12.5	15	15



Thorax

CM	2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
kVp	68	72	76	80	84	88	92	96	100	100	100	100	100	100	100
mAs	5	5	5	5	5	5	5	7.5	7.5	7.5	7.5	10	10	12.5	12.5



Spine, Skull, Pelvis

CM	2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
kVp	62	64	66	68	70	74	76	80	84	90	94	100	100	100	100
mAs	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	12.5	15	15



Extremities/Exotics

CM	2	4	6	8	10
kVp	72	72	72	76	76
mAs	5	5	5	5	5

Note: If images is too "grainy" INCREASE mAs 50%
If image is too "saturated" DECREASE mAs 50%



800.819.5538 soundvet.com
resources.soundvet.com