

SmartDR Premier Technique Guideline

Weight (lbs)	Extremity/Skull	Other
1 - 20 lbs	65 kVp / 4.0 mAs	65 kVp / 4.0 mAs
21 - 45 lbs	65 kVp / 4.0 mAs	75 kVp / 4.0 mAs
46 - 80 lbs	65 kVp / 4.0 mAs	90 kVp / 5.0 mAs
81 - 100 lbs	65 kVp / 4.0 mAs	95 kVp / 6.0 mAs
100+ lbs	90 kVp / 6.0 mAs	100 kVp / 6.5 mAs