

## Sound SmartDR+ CsI Technique Chart

	Thickness (cm)	1-10	11-15	16-25	26-30+
	Skull, Spine, Pelvis, Extremities	kVp	66	76	86
	mAs	2.5	2.5	2.5	5.0
Thorax Lateral	kVp	70	78	80	90
	mAs	1.2	1.2	1.7	2.5
Thorax VD	kVp	66	76	86	96
	mAs	1.5	2.0	2.7	2.7
Abdomen Lateral	kVp	65	70	78	84
	mAs	2.0	2.6	2.6	2.6
Abdomen VD	kVp	80	86	92	98
	mAs	2.0	2.6	3.0	4.6





